

Art Unit: ***

Claims PTO

12/15/2004

AMW

BEST AVAILABLE COPY

1. A device comprising a sliding member and a fixed member, the sliding member indicative of the moisture condition of a wound and the fixed member indicative of the healing phase condition of a wound and containing instructions as to the wound treatment regiment resulting from the alignment of the sliding member's particular wound moisture condition with the fixed member's particular wound healing phase condition.

2. (Amended) A device comprising a sliding member and a fixed member, wherein the sliding member is indicative of the wound's healing phase condition and the fixed member is indicative of the wound's moisture condition.

3. A device useful for prescribing a treatment regiment for chronic wounds comprising:

(a) a base comprising first and second elongated fixed members, the members being spaced to receive a moveable slide;

(b) the moveable slide being located between the two spaced elongated fixed members;

(c) the first of the elongated fixed members separated into regions indicating the healing phase

Art Unit: ***

condition of the wound in terms of whether the wound is in the necrotic condition; fibrinous slough or infection condition; granulation condition; or epithelialization condition; each of these healing phase condition regions further containing indicators corresponding to the moisture condition of the wound indicating whether the moisture level of the wound is in the wet condition; moist condition, or dry condition;

(d) the moveable slide being separated into three regions corresponding to the moisture condition of the wound in terms of whether the wound is in the wet condition, moist condition, or dry condition; each of these moisture condition regions further containing wound healing phase indicators corresponding to the healing phase condition of the wound in terms of whether the wound is in the necrotic condition; fibrinous slough or infection condition, granulation condition; or epithelialization condition; and

(e) the second of the elongated fixed members comprising treatment descriptors comprising instructions for wound treatment regiments corresponding to the wound phase healing and moisture conditions which result from the alignment of the healing and moisture condition indicators of the first elongated fixed member with the healing phase and moisture condition indicators of the moveable slide.

Art Unit: ***

4. The device of claim 3, wherein the wound healing phase condition regions are of the first elongated fixed member shaded or colored to distinguish the regions from among themselves.

5. The device of claim 4, wherein the wound moisture condition regions of the movable slide are shaded or colored to distinguish the regions from among themselves.

6. The device of claim 5, wherein the healing phase condition regions have the following coloring scheme of black to depict the necrotic condition, yellow to depict the fibrinous slough or infection condition, red to depict the granulation condition and pink to depict the epithialization condition.

7. The device of claim 6, wherein the wound moisture condition regions have the following coloring scheme of blue to depict a wet wound, green to depict a moist wound, and yellow to depict a dry wound.

8. The device of claim 7, wherein the moisture condition indicators of the first elongated fixed member, the wound healing phase indicators of the moveable slide, and the indicators of the second elongated fixed member are shaded or colored.

Art Unit: ***

9. The device of claim 8, wherein the moisture condition indicators of the first elongated member, the wound healing phase indicators of the sliding member, and the indicators of the second elongated member each individually comprise two side-by-side colored bars.

10. The device of claim 9, wherein the indicators comprise the following left to right color scheme for the indicated wound healing and moisture conditions:

(a) FIRST ELONGATED FIXED MEMBER:

(i) Necrotic Condition:	<u>Left Bar</u>	<u>Right Bar</u>
wet (1101)	black	blue
moist (1102)	black	green
dry (1103)	black	yellow
(ii) Fibrinous Slough/Infection:	<u>Left Bar</u>	<u>Right Bar</u>
wet (1201)	yellow	blue
moist (1202)	yellow	green
dry (1203)	yellow	yellow
(iii) Granulation Condition:	<u>Left Bar</u>	<u>Right Bar</u>
wet (1301)	red	blue
moist (1302)	red	green
dry (1303)	red	yellow

Art Unit: ***

(iv) Epithelialization Condition:	<u>Left Bar</u>	<u>Right Bar</u>
wet (1401)	pink	blue
moist (1402)	pink	green
dry (1403)	pink	yellow;
(b) MOVEABLE SLIDE:		
(i) Wet Condition:	<u>Left Bar</u>	<u>Right Bar</u>
necrotic (2101)	blue	black
fibrinous slough/infection (2102)	blue	yellow
granulation (2103)	blue	red
epithelialization (2104)	blue	pink
(ii) Moist Condition:	<u>Left Bar</u>	<u>Right Bar</u>
necrotic (2201)	green	black
fibrinous slough/infection (2202)	green	yellow
granulation (2203)	green	red
epithelialization (2204)	green	pink
(iii) Dry Condition:	<u>Left Bar</u>	<u>Right Bar</u>
necrotic (2201)	yellow	black
fibrinous slough/infection (2202)	yellow	yellow
granulation (2203)	yellow	red
epithelializatoin (2204)	yellow	pink; and

Art Unit: ***

(c) SECOND ELONGATED FIXED MEMBER:

(i) Necrotic Condition:	<u>Left Bar</u>	<u>Right Bar</u>
wet (3101)	black	blue
moist (3103)	black	green
dry (3103)	black	yellow
(ii) Fibrinous Slough/Infection:	<u>Left Bar</u>	<u>Right Bar</u>
wet (3201)	yellow	blue
moist (3202)	yellow	green
dry (3203)	yellow	yellow
(iii) Granulation Condition:	<u>Left Bar</u>	<u>Right Bar</u>
wet (3301)	red	blue
moist (3302)	red	green
dry (3303)	red	yellow
(iv) Epithelialization Condition:	<u>Left Bar</u>	<u>Right Bar</u>
wet (3401)	pink	blue
moist (3402)	pink	green
dry (3403)	pink	yellow.

11. A method for prescribing a treatment regiment
for a wound comprising the steps of:

(a) assessing the wound healing phase condition;

(b) assessing the wound moisture condition;

Art Unit: ***

(c) providing a device comprising at least two fixed members and a sliding member, the first of the fixed members representing the wound healing condition, the sliding member representing the wound moisture condition, and the second of the fixed members comprising wound treatment instructions;

(d) aligning the moisture condition of the sliding member with the wound healing phase condition of the first fixed member corresponding to the assessed wound healing phase and moisture condition and then prescribing the wound treatment regimen indicated on the second fixed member resulting from the alignment of the sliding member and the first fixed member.